DCM025-12
Ed. 01/2015

PROGESTERONE SALIVA ELISA

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Ed. 01/2015

DCM025-12

DiaMetra S.r.l. Headquater: Via Calabria 15,
20090 SEGRATE (MI) Italy
Tel. +39-02-2139184
Fax +39-02-2133354

Manufactory: Via Pozzuolo 14,
06038 SPELLO (PG) Italy
Tel. +39-0742-24851
Fax +39-0742-316197
E-mail: info@diametra.com

DCM025-12
Ed. 01/2015

PROGESTERONE SALIVA ELISA

for routine analysis


Direct immunoenzymatic determination of Progesterone in saliva.

IVD

LOT

See external label

Σ = 96 tests

REF DKO025

INTENDED USE

Competitive immunoenzymatic colorimetric method for quantitative determination of Progesterone concentration in saliva.

Progesterone Saliva ELISA kit is intended for laboratory use only.

1. CLINICAL SIGNIFICANCE

Progesterone is a C-21 steroid hormone involved in the female menstrual cycle, pregnancy (supports *gestation*) and embryogenesis of humans and other species. Progesterone is the major naturally occurring human progestagen. Progesterone is important for aldosterone (mineralocorticoid) synthesis, as 17-hydroxyprogesterone is for cortisol (glucocorticoid).

Progesterone levels are relatively low in children and postmenopausal women. Adult males have levels similar to those in women during the follicular phase of the menstrual cycle. In women, progesterone levels are relatively low during the preovulatory phase of the menstrual cycle, rise after ovulation, and are elevated during the luteal phase. If pregnancy occurs, progesterone levels are maintained at luteal levels initially. After delivery of the placenta and during lactation, progesterone levels are very low. The fall in progesterone levels following delivery is one of the triggers for milk production. Progesterone is produced in the adrenal glands, the gonads (specifically after ovulation in the corpus luteum), the brain, and, during pregnancy, in the placenta. Progesterone converts the endometrium to its secretory stage to prepare the uterus for implantation. If pregnancy does not occur, progesterone levels will decrease, leading, in the human, to menstruation. Progesterone belongs to the group of neurosteroids that are found in high concentrations in certain areas in the brain and are synthesized there. Neurosteroids affect synaptic functioning, are neuroprotective, and affect myelination. Progesterone has multiple effects outside of the reproductive system. Progesterone is thermogenic, it reduces spasm and relaxes smooth muscle. Bronchi are widened and mucus regulated. Progesterone acts as an antiinflammatory agent and regulates the immune response. Progesterone also assists in thyroid function, in bone building by osteoblasts. Measurement of serum progesterone concentrations have been used in evaluating ovarian function.

2. PRINCIPLE

Progesterone (antigen) in the sample competes with progesterone conjugated with horseradish peroxidase (enzyme-labelled antigen) for binding to the limited number of anti-progesterone antibody sites on the microplates (solid phase).

After incubation, the bound/free separation is performed by a simple solid-phase washing.

The enzyme in the bound-fraction reacts with the Substrate (H_2O_2) and the TMB Substrate and develops a blu color that changes into yellow when the Stop Solution (H_2SO_4) is added.

The colour intensity is inversely proportional to the Progesterone concentration in the sample.

Progesterone concentration in the sample is calculated through a calibration curve.

3. REAGENT, MATERIAL AND INSTRUMENTATION

3.1. Reagent and material supplied in the kit

1. Calibrators (7 vials, 1 mL each)

CAL0	REF DCE002/2506-0
CAL1	REF DCE002/2507-0
CAL2	REF DCE002/2508-0
CAL3	REF DCE002/2509-0
CAL4	REF DCE002/2510-0
CAL5	REF DCE002/2511-0
CAL6	REF DCE002/2512-0

2. Controls (2 vials, 1 mL each)

Low Control	REF DCE045/2504-0
High Control	REF DCE045/2505-0

3. Conjugate (1 vial, 26 mL)

Progesterone conjugated with horseradish peroxidase (HRP) REF DCE002/2502-0

4. Coated Microplate (1 microplate breakable)

Policlonal anti Progesterone antibody adsorbed on microplate REF DCE002/2503-0

5. TMB Substrate (1 vial, 25 mL)

H_2O_2 -TMB (0,26 g/L) (avoid any skin contact) REF DCE004/2504-0

6. Stop Solution (1 vial, 14 mL)

Sulphuric acid 0.5 mol/L (avoid any skin contact) REF DCE005/2505-0

7. 40X Conc. Wash Solution (1 vial, 30 mL)

Buffered solution REF DCE057/2557-0

3.2. Reagents necessary not supplied

Distilled water.

3.3. Auxiliary materials and instrumentation

Automatic dispenser.

Microplates reader (450 nm, 620-630 nm)

Saliva Collection Device REF DKO063

Note

Store all reagents at 20^o8°C in the dark.

Open the bag of reagent 4 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened the microplate is stable until the expiry date of the kit.

Do not remove the adhesive sheets on the unutilised strips.

4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the reagent should be handled in the same manner as potentially infectious material.
- Some reagents contain small amounts of Proclin 300^R as preservative. Avoid the contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to directed sunlight, metals or oxidants. Do not freeze the solution.
- This method allows the determination of Progesterone from 10 pg/mL to 5000 ng/mL.
- The clinical significance of the Progesterone determination can be invalidated if the patient was treated with cortisone or natural or synthetic steroids.

5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated. Opened kits retain activity for two months if stored as described above.

- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled saliva samples.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Plate readers measure vertically. Do not touch the bottom of the wells.

6. PROCEDURE

6.1. Preparation of the Calibrators (C₀...C₆)

The Calibrators are ready to use and have the following concentration:

	C ₀	C ₁	C ₂	C ₃	C ₄	C ₅	C ₆
pg/mL	0	10	50	100	500	1000	5000

For SI Units: pg/mL x 3.18 = pmol/L

6.2. Preparation of Wash Solution

Dilute the content of each vial of the "40X Wash Solution" with distilled water to a final volume of 1200 mL prior to use. For smaller volumes respect the 1:40 dilution ratio. The diluted wash solution is stable for 2 weeks at room temperature.

6.3. Preparation of Conjugate

The Conjugate is ready for use.

6.4. Preparation of the Sample

This kit allows determination of Progesterone in saliva samples.

It is recommended to collect saliva samples with a centrifuge glass tube and a plastic straw, or with the Diametra *Saliva Collection Device* (DKO063). Other commercially available sample collector devices have not been tested.

6.4.1. Method and Limitations

Collect saliva samples at the times indicated.

If no specific instructions have been given the saliva samples may be collected at any time; for saliva collection, the following should be noted:

- a) If saliva collection is carried out in the morning ensure that this is carried out prior to brushing teeth
- b) During the day allow 1 hour after a meal, oral intake of pharmaceutical drugs or tooth cleaning.
- c) It is very important that a good clear sample is received – i.e. no contamination with food, lipstick, blood (bleeding gums) or other extraneous materials.

6.4.2. Saliva Processing Instructions with Saliva Collection Device Diametra

- 1) Let the saliva flow down through the straw into the centrifuge glass tube.
- 2) Centrifuge the sample for 15 minutes at 3000 rpm
- 3) Store at – 20°C for at least 1 hour
- 4) Centrifuge again for 15 minutes at 3000 rpm
- 5) The saliva sample is now ready to be tested.
- 6) Store the sample at 2-8°C for one week or at – 20°C for longer time.

6.5. Procedure

- **Allow all reagents to reach room temperature (22-28°C).** At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₆), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample	Blank
Sample		100 µL	
Calibrator C ₀ -C ₅	100 µL		
Conjugate	200 µL	200 µL	
Gently mix the plate for 10 seconds and incubate at room temperature (22±28°C) for 1 hour. Remove the content from each well; wash the wells 3 times with 400 µL of diluted wash solution. Strike the inverted well sharply on absorbent paper towel to remove residual droplets. Important note: the sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!			
TMB Substrate	200 µL	200 µL	200 µL
Incubate at room temperature (22±28°C) for 15 minutes in the dark			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

7. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of Progesterone for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the calibration curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

8. RESULTS

8.1. Calculation of Results

1. Calculate the average absorbance values for each set of Calibrators, Controls and patient samples.
2. Using linear graph paper, construct a calibration curve by plotting the mean absorbance obtained from each Calibrator against its concentration with absorbance value on the vertical (Y) axis and concentration on the horizontal (X) axis.
3. Using the mean absorbance value for each sample determine the corresponding concentration from the calibration curve.
4. Automated method: the results in the IFU have been calculated automatically using a 4 PL (4 Parameter Logistics) curve fit. 4 Parameter Logistics is the preferred calculation method. Other data reduction functions may give slightly different results.
5. The concentration of the samples can be read directly from this calibration curve. Samples with concentrations higher than that of the highest Calibrator have to be further diluted or reported as > 5000 pg/mL. For the calculation of the concentrations this dilution factor has to be taken into account.

8.2. Example of a typical calibration curve

The following data is for demonstration only and cannot be used in place of data generation at the time of assay.

Calibrator	Absorbance Units (450 nm)
Calibrator 0	1.890
Calibrator 1	1.710
Calibrator 2	1.540
Calibrator 3	1.390
Calibrator 4	0.950
Calibrator 5	0.710
Calibrator 6	0.440

9. EXPECTED NORMAL VALUES

In order to determine the normal range of Progesterone in saliva, saliva samples from 80 adult male and 120 female apparently healthy subjects, ages 21 to 75 years, were collected in the morning and analyzed using the Diametra Progesterone ELISA kit.

	Age group	Progesterone saliva (pg/mL)
Women	21 - 50 years; n = 40 Follicular phase n = 40	19,6 - 86,5
	21 - 50 years; n = 40 Luteal phase	99,1 - 332,6
	21 - 50 years; n = 40 Postmenopausal	6,0 - 56,4
Men	21 - 50 years; n = 40	12,7 - 57,4
	41 - 75 years; n = 40	15,2 - 65,1

Please pay attention to the fact that the determination of a range of expected values for a "normal" population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

10. PERFORMANCE AND CHARACTERISTICS

10.1. Precision

10.1.1. Intra Assay Variation

The intra-assay variation was determined by 20 replicate measurements of 5 saliva samples within one run.

Mean (pg/mL)	1328.7	650.1	293.8	186.9	23.3
SD	75.9	30.4	15.8	10.6	1.8
CV %	5.7	4.7	5.4	5.7	7.6
n =	20	20	20	20	20

10.1.2. Inter Assay Variation

The inter-assay (between-day) variation was determined by duplicate measurements of 5 saliva samples over 10 days.

Mean (pg/mL)	60.4	630.3	1574.8	360.3	635.1
SD	4.3	37.6	64.1	28.7	33.4
CV %	7.1	6.0	4.1	8.0	5.3
n =	9	9	9	9	9

10.2. Recovery

Recovery of the Diametra Progesterone ELISA was determined by adding increasing amounts of the analyte to 5 different saliva samples containing different amounts of endogenous analyte. Each sample (native and spiked) was assayed and analyte concentrations of the samples were calculated from the calibration curve. The percentage recoveries were determined by comparing expected and measured values of the samples.

	1	2	3	4	5
Conc.	13.9	87.9	147.8	283.0	386.8
% Recovery	96.4	103.2	100.3	97.3	102.2
Range Recovery	91.1 101.6	101.1 104.8	94.9 103.8	92.7 106.6	97.8 105.7

10.3. Linearity

In total six saliva samples containing different amounts of analyte were serially diluted with the Calibrator 0 and assayed with the Diametra Progesterone saliva ELISA. Three of these samples were serially diluted directly, and the other 3 samples at first were spiked with Progesterone and the serially diluted up to 1:128. The percentage recovery was calculated by comparing the expected and measured values for Progesterone. A linearity of 3,8 - 4600 pg/mL has been identified as

the usable range for this assay. Samples above this range must be diluted and re-run.

	1	2	3	4	5	6
Conc.	58.3	98.7	1073	2802	6000	5000
% Recovery	104.1	97.5	95.79	98.96	102.6	103.5
Range Recovery	99.9	89.9	91.2	93.6	99.1	94.4
	106.7	100.8	105.5	107.3	106.5	107.7

10.4. Sensitivity

The lowest detectable concentration of Progesterone that can be distinguished from the Calibrator 0 is 3.8 pg/mL at the 95% confidence limit.

10.5. Specificity

The following materials have been evaluated for cross reactivity. The percentage indicates cross reactivity at 50% displacement compared to Progesterone.

Steroid	% cross reactivity
Progesterone	100
Desoxycorticosterone	1.1
Pregnenolone	0.35
17a-Hydroxyprogesterone	0.3
Corticosterone	0.2
11-Desoxycortisol	0.1
Estriol	< 0.1
17b Estradiol	< 0.1
Testosterone	< 0.1
Cortisone	< 0.1
DHEA-S	< 0.02
Cortisol	< 0.02

10.6. Comparison studies

Diametra Progesterone saliva ELISA kit was compared to a commercially available LIA Progesterone saliva assay. 306 saliva samples were collected from adult men and women ages 20-75.

The linear regression curve was calculated:

$$y = 0,912x + 6,066$$

$$r^2 = 0,937$$

An additional study was performed using 101 saliva samples from adult men and women ages 20-75. These samples were compared to the LIA method and the linear regression curve was calculated:

$$y = 0,934x + 2,233$$

$$r^2 = 0,992$$

10.7. Interfering substances

- Blood contamination of more than 0,16% in saliva samples will affect results, and usually can be seen by eye.
- Concentrations of Sodium Azide 0,02% interferes in this assay and may lead to false results.

11. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

BIBLIOGRAPHY

- Wisdom, G.B. *Clin. Chem.* 22 (8), 1243- 1255 (1976)
- De Villa, G.O., et al *J. Clin. Endocr. Metab.* 35, 458 (1972)
- Joyce, B.G., et al *Steroids* 29, no 6, 761, (1977)
- Winkel P., et al *Clin. Chem.* 22 (4), 422 1976
- Rajkowski K. N., et al *Steroids* 29, no 5 (1977)
- Y. Lu, et al *Fertility and Sterility* Vol.71 N.5 863 (1999)
- Miller W.L. (1988): Molecular biology of steroid hormone synthesis, *Endocrin. Rev.* 9, 295-318
- Filicori M., Butler J.P., Crowley W.F. Jr. (1984): Neuroendocrine regulation of the corpus luteum in the human, *J.Clin. Invest.* 73, 1638
- Csapol Al., Pulkkinen M. O., Wiest W.G. (1973): Effects of lutectomy and progesterone replacement therapy in early pregnancy patients. *Am J. Obstet. Gynecol.* 115, 759
- Henson M.C. (1998): Pregnancy maintenance and the regulation of placental progesterone biosynthesis in the baboon. *Human Reproduction Update*, 4, 389-405

Ed. 01/2015

DCM025-12

DiaMetra S.r.l. Headquarter: Via Calabria 15,
20090 SEGRATE (MI) Italy

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


ERROR POSSIBLE CAUSES / SUGGESTIONS

No colorimetric reaction

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

Too low reaction (too low ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

Too high reaction (too high ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

Unexplainable outliers

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation